

ACR Ambulance Child Restraint



ParAid's ACR (Ambulance Child Restraint) is an innovative, flexible and fully adjustable harnessing system for the safe and effective transport of infants and children in an ambulance, and is colour coded for easy selection in 3 sizes from 5kg to 45kg.

Quick release clips dock with the ACR harness, holding the patient in place to prevent potentially dangerous movement during transportation. This product is manufactured under an exclusive licence agreement.

- BS EN 1865 & 1789 compliant
- Weight range 5kg - 45kg
- Universally compatible with any stretcher
- Smart colour coding for easy selection of 3 sizes
- Machine washable
- Compact storage bag
- Quick release chest strap for medical intervention
- Breathable fabric maximises airflow and meets infection control standards
- Allows rapid transition from sitting to flat in an emergency or to recovery position
- Padded for patient comfort and reduced risk of pressure sores during long transfers
- Can be used with stretcher back rest in raised position



ParAid House, Weston Lane, Birmingham B11 3RS T: +44 (0)121 700 7455 F: +44 (0)121 700 7454

E: info@paraid.co.uk W: www.paraid.com